STANDARD Q

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
STANDARD™ Q COVID-19 Ag Test

SD BIOSENSOR

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

\* It is recommended to perform ANTIGEN TEST WITHIN 9DAYS after symptom onset.


## KIT CONTENTS

## PREPARATION AND TEST PROCEDURE

## [Preparation]

- Carefully read instructions for using the STANDARD Q COVID-19 Ag Test.
- Check the expiry date at the back of the foil pouch.  
Do not use the kit, if expiry date has passed.
- Check the test device and the desiccant pack in the foil pouch.

## [Test Procedure]

## 1 Nasopharynx method


- Insert a sterile swab into the nostril of the patient, reaching the surface of the posterior nasopharynx.

- Swab over the surface of the posterior nasopharynx for 5 seconds.


- Withdraw the sterile swab from the nasal cavity.


## 2 Nasal cavity method


- Insert a sterile swab about 2cm along the patient's nasal cavity.
- Rotate the sterile swab about 5 times, then withdraw the sterile swab from the nasal cavity.

- Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.

- Apply 3 drops of extracted specimen to the specimen well of the test device.
- Read the test result in 15-30 minutes.

CAUTION

- Do not read test results after 30 minutes. It may give false results.

## INTERPRETATION OF TEST RESULT

\* "C" Control Line | "T" Test Line

Negative	
Positive	
Invalid	

- A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).
  - A colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T).
  - Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.
- \* The presence of any line no matter how faint the result is considered positive.
- \* Positive results should be considered in conjunction with the clinical history and other data available.

## KIT STORAGE AND STABILITY

Store the kit at 2-30° C / 36-86° F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

## WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the extraction buffer tube of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.