

Nitrile Gloves

NITRILE GLOVES POWDER FREE



PHYSICAL DIMENSIONS

| DIMENSIONS | Standards | |
|--------------------------------|--|---|
| | VRG KHAI HOAN | ASTM D6319 |
| Length (mm) | 230 min | 220 min (XS, S) 230 min (M, L, XL) |
| Width (mm) | 75 ± 5 (XS) 85 ± 5 (S) 95 ± 5 (M) 105 ± 5 (L) 115 ± 5 (XL) | 70 ± 10 (XS) 80 ± 10 (S) 95 ± 10 (M) 110 ± 10 (L) 120 ± 10 (XL) |
| Thickness- Single wall (mm) | Fingers : 0.08 mm min Palm : 0.06 mm min | Fingers : 0.050 mm min Palm : 0.05 mm min |

PHYSICAL PROPERTIES AND BIOCOMPATIBILITY

| | | |
|-----------------|--|--|
| Tensile | <p>Tensile strength (MPA) Before aging: 18Mpa min After aging: 20Mpa min</p> <p>Elongation at break (%) Before aging: 600% min After aging: 500% min</p> | <p>Tensile strength (MPA) Before aging: 14Mpa min After aging: 14Mpa min</p> <p>Elongation at break (%) Before aging: 500% min After aging: 400% min</p> |
| Powder Content | 2 mg/glove maximum | |
| Protein Content | Free Protein | |

